



2022

Product Catalogue



Idea for Diagnosis

BiON
IDEA FOR LAB

(주)바이온라이프사이언스

E. bion@bionlifescience.com
www.bionlifescience.com

Bion Lifescience, Inc. registered to various certification



Elastomer Type

For the collection of clinical specimens containing virus, germs



BION™ VIROSWAB(Sterile) and BION™ nVIROSWAB(Non-Sterile) features

- Head material is Elastomer.
- There are two head types.(Flower type, Screw type)
- BION™ VIROSWAB and BION™ nVIROSWAB are manufactured in ISO 13485 certified facility.
- Medical devices registered with the Korea MFDS
- Applied for a patent.

Elastomer Swab and Nylon Flocked Swab Comparison

Criteria	Type	Elastomer Swab	Nylon Flocked Swab
Head's Material		Elastomer	Nylon Fiber
Bonding		Non	water soluble/ fat-soluble bonding

BION™ VIROSWAB(Sterile) Specification

	ENS-F81-S	ENS-F91-S	ENS-S81-S	ENS-S91-S	EOS-F79-S	EOS-F91-S	EOS-S79-S	EOS-S91-S	ENOS-F8179-S	ENOS-F9191-S	ENOS-F8179-S	ENOS-S9191-S
Use	Nasal				Oral				Nasal & Oral			
Head Type	Flower		Screw		Flower		Screw		Flower		Screw	
Over-all length	153mm											
Break Point	81mm	91mm	81mm	91mm	79mm	91mm	79mm	91mm	81mm & 79mm	91mm & 91mm	81mm & 79mm	91mm & 91mm


BION™ nVIROSWAB(Non-Sterile) Specification

	ENS-F81	ENS-F91	ENS-S81	ENS-S91	EOS-F79	EOS-F91	EOS-S79	EOS-S91	ENOS-F8179	ENOS-F9191	ENOS-F8179	ENOS-S9191
Use	Nasal				Oral				Nasal & Oral			
Head Type	Flower		Screw		Flower		Screw		Flower		Screw	
Over-all length	153mm											
Break Point	81mm	91mm	81mm	91mm	79mm	91mm	79mm	91mm	81mm & 79mm	91mm & 91mm	81mm & 79mm	91mm & 91mm

BION™ VIROSWAB(Sterile) and BION™ nVIROSWAB(Non-Sterile) applications

- Break point 81mm and 79mm swabs can be used with 5ML Screw Tubes and RDT (POCT).
- Break point 91mm swabs can be used with 12ML Screw Tubes.

Instruction for use



Tilt the patient's head a little backward.

1-1. Oral swab

Without touching the sides of the mouth, rub the back of throat and the tonsil area on both side of the throat.

1-2. Nasal Swab

Insert the tip horizontally until it reaches to the Nasopharynx and rotate the swab 7 times.

2. Take out the swabs carefully.

3. Insert the swabs into a container until it reaches to the bottom.

4. Break the swabs at the molded break point.

5. Close the cap firmly to ensure no leakage.

BION™ VIROSWAB

Specimen Collection Swab
Elastomer Type

Nasal Type

Head Part



Flower Type



Screw Type

ENS-F81-S



ENS-F91-S



ENS-S81-S



ENS-S91-S



Model	Use	Head Type	Thickness of Head	Length of Head	Length of Head to Break Point	Over-all length	Sterile
ENS-F81-S	Nasal	Flower	4mm	20mm	81mm	153mm	0
ENS-F91-S	Nasal	Flower	4mm	20mm	91mm	153mm	0
ENS-S81-S	Nasal	Screw	4mm	20mm	81mm	153mm	0
ENS-S91-S	Nasal	Screw	4mm	20mm	91mm	153mm	0

Specimen Collection Swab
Elastomer Type

Oral Type

Head Part



Flower Type



Screw Type

EOS-F79-S



EOS-F91-S



EOS-S79-S



EOS-S91-S



Model	Use	Head Type	Thickness of Head	Length of Head	Length of Head to Break Point	Over-all length	Sterile
EOS-F79-S	Oral	Flower	5mm	20mm	79mm	153mm	0
EOS-F91-S	Oral	Flower	5mm	20mm	91mm	153mm	0
EOS-S79-S	Oral	Screw	5mm	20mm	79mm	153mm	0
EOS-S91-S	Oral	Screw	5mm	20mm	91mm	153mm	0


Specimen Collection Swab Elastomer Type

Nasal & Oral Type

Head Part



Flower-Nasal Type



Flower-Oral Type

ENOS-F8179-S



ENOS-F9191-S



Model	Use	Head Type	Thickness of Head	Length of Head	Length of Head to Break Point	Over-all length	Sterile
ENOS-F8179-S	Nasal & Oral	Flower	Nasal: 4mm / Oral: 5mm	20mm / 20mm	Nasal : 81mm / Oral : 79mm	153mm / 153mm	0
ENOS-F9191-S	Nasal & Oral	Flower	Nasal: 4mm / Oral: 5mm	20mm / 20mm	Nasal : 91mm / Oral : 91mm	153mm / 153mm	0

Head Part



Screw-Nasal Type



Screw-Oral Type

ENOS-S8179-S



ENOS-S9191-S



Model	Use	Head Type	Thickness of Head	Length of Head	Length of Head to Break Point	Over-all length	Sterile
ENOS-S8179-S	Nasal & Oral	Screw	Nasal: 4mm / Oral: 5mm	20mm / 20mm	Nasal : 81mm / Oral : 79mm	153mm / 153mm	0
ENOS-S9191-S	Nasal & Oral	Screw	Nasal: 4mm / Oral: 5mm	20mm / 20mm	Nasal : 91mm / Oral : 91mm	153mm / 153mm	0

BION™ nVIROSWAB

Specimen Collection Swab Elastomer Type

Nasal Type

Head Part



Flower Type



Screw Type

ENS-F81



ENS-F91



ENS-S81



ENS-S91



Model	Use	Head Type	Thickness of Head	Length of Head	Length of Head to Break Point	Over-all length	Sterile
ENS-F81	Nasal	Flower	4mm	20mm	81mm	153mm	X
ENS-F91	Nasal	Flower	4mm	20mm	91mm	153mm	X
ENS-S81	Nasal	Screw	4mm	20mm	81mm	153mm	X
ENS-S91	Nasal	Screw	4mm	20mm	91mm	153mm	X

Specimen Collection Swab Elastomer Type

Oral Type

Head Part



Flower Type



Screw Type

EOS-F79



EOS-F91



EOS-S79



EOS-S91



Model	Use	Head Type	Thickness of Head	Length of Head	Length of Head to Break Point	Over-all length	Sterile
EOS-F79	Oral	Flower	5mm	20mm	79mm	153mm	X
EOS-F91	Oral	Flower	5mm	20mm	91mm	153mm	X
EOS-S79	Oral	Screw	5mm	20mm	79mm	153mm	X
EOS-S91	Oral	Screw	5mm	20mm	91mm	153mm	X


Specimen Collection Swab Elastomer Type

Nasal & Oral Type

Head Part



Flower-Nasal Type



Flower-Oral Type

ENOS-F8179



ENOS-F9191



Model	Use	Head Type	Thickness of Head	Length of Head	Length of Head to Break Point	Over-all length	Sterile
ENOS-F8179	Nasal & Oral	Flower	Nasal: 4mm / Oral: 5mm	20mm / 20mm	Nasal : 81mm / Oral : 79mm	153mm / 153mm	X
ENOS-F9191	Nasal & Oral	Flower	Nasal: 4mm / Oral: 5mm	20mm / 20mm	Nasal : 91mm / Oral : 91mm	153mm / 153mm	X

Head Part



Screw-Nasal Type



Screw-Oral Type

ENOS-S8179



ENOS-S9191



Model	Use	Head Type	Thickness of Head	Length of Head	Length of Head to Break Point	Over-all length	Sterile
ENOS-S8179	Nasal & Oral	Screw	Nasal: 4mm / Oral: 5mm	20mm / 20mm	Nasal : 81mm / Oral : 79mm	153mm / 153mm	X
ENOS-S9191	Nasal & Oral	Screw	Nasal: 4mm / Oral: 5mm	20mm / 20mm	Nasal : 91mm / Oral : 91mm	153mm / 153mm	X

Screw Caps and Tubes 2ML

For the collection of clinical specimens containing virus, germs And for holding reagent



Natural Cap
& Clear Tube



Amber Cap
& Amber Tube

Our Screw Caps and Tubes make your choice best.

- Good design for easy grab.
- Self Standing makes your hand free.
- Neat and clear appearance.
- Various color choice.

Screw Caps and Tubes 2ML

For the collection of clinical specimens containing virus, germs And for holding reagent



O-ring is not included with this product

Our cap design functions well without O-ring.

- Our Screw caps provide a tight seal.
 - There is no 'O-ring' so there is no risk of contamination from rubber particles.
 - Easy $\frac{3}{4}$ turn to open and close, and It turns smoothly.
- Our caps give you convenience.

Screw Caps and Tubes 2ML

For the collection of clinical specimens containing virus, germs And for holding reagent

2ML Screw Caps and Tubes, Non-Sterile

Cat. No.	Description	Cap Colors	Tube Color	Dimensions Ø x h	Packaging
SCT-2ML-SS-RC	Screw Caps and Tubes Non-Sterile Conical bottom, Self Standing	Red	Clear	13mm x 47mm	6000 Tubes & Caps (1000 Tubes x 6pks & 3000 Caps x 2pks)
SCT-2ML-SS-OC		Orange	Clear		
SCT-2ML-SS-YC		Yellow	Clear		
SCT-2ML-SS-GC		Green	Clear		
SCT-2ML-SS-BC		Blue	Clear		
SCT-2ML-SS-PUC		Purple	Clear		
SCT-2ML-SS-WC		White	Clear		
SCT-2ML-SS-AMC		Amber	Clear		
SCT-2ML-SS-NAC		Natural	Clear		
SCT-2ML-SS-AMX		Amber	Amber		



2ML Screw Caps and Tubes, Sterile

Cat. No.	Description	Cap Colors	Tube Color	Dimensions Ø x h	Packaging
SCT-2ML-SS-RC-S	Screw Caps and Tubes Sterile by E-beam Irradiation Conical bottom, Self Standing	Red	Clear	13mm x 47mm	6000 Tubes & Caps (1000 Tubes x 6pks & 3000 Caps x 2pks)
SCT-2ML-SS-OC-S		Orange	Clear		
SCT-2ML-SS-YC-S		Yellow	Clear		
SCT-2ML-SS-GC-S		Green	Clear		
SCT-2ML-SS-BC-S		Blue	Clear		
SCT-2ML-SS-PUC-S		Purple	Clear		
SCT-2ML-SS-WC-S		White	Clear		
SCT-2ML-SS-AMC-S		Amber	Clear		
SCT-2ML-SS-NAC-S		Natural	Clear		
SCT-2ML-SS-AMX-S		Amber	Amber		




Screw Caps and Tubes 2ML

For the collection of clinical specimens containing virus, germs And for holding reagent

2ML Screw Caps, Non-Sterile

Cat. No.	Description	Colors	Dimensions Ø x h	Packaging
SC-2ML-R	Screw Caps Non-Sterile 	Red	13mm x 7.4mm	Custom order
SC-2ML-O		Orange		
SC-2ML-Y		Yellow		
SC-2ML-G		Green		
SC-2ML-B		Blue		
SC-2ML-PU		Purple		
SC-2ML-W		White		
SC-2ML-AM		Amber		
SC-2ML-NA		Natural		

2ML Screw Caps, Sterile

Cat. No.	Description	Colors	Dimensions Ø x h	Packaging
SC-2ML-R-S	Screw Caps Sterile by E-beam Irradiation 	Red	13mm x 7.4mm	Custom order
SC-2ML-O-S		Orange		
SC-2ML-Y-S		Yellow		
SC-2ML-G-S		Green		
SC-2ML-B-S		Blue		
SC-2ML-PU-S		Purple		
SC-2ML-W-S		White		
SC-2ML-AM-S		Amber		
SC-2ML-NA-S		Natural		

Screw Caps and Tubes 2ML

For the collection of clinical specimens containing virus, germs And for holding reagent

2ML Screw Tubes

Cat. No.	Description	Colors	Dimensions Ø x h	Packaging
ST-2ML-SS-C	Conical bottom, Self Standing, Non-Sterile	Clear	12.5mm x 45.5mm	Custom order
ST-2ML-SS-C-S	Conical bottom, Self Standing, Sterile	Clear		
ST-2ML-SS-AM	Conical bottom, Self Standing, Non-Sterile	Amber		
ST-2ML-SS-AM-S	Conical bottom, Self Standing, Sterile	Amber		



Clear



Amber

Color List



Red



Orange



Yellow



Green



Natural

Non-Sterile	SC-2ML-R	SC-2ML-O	SC-2ML-Y	SC-2ML-G	SC-2ML-NA
Sterile	SC-2ML-R-S	SC-2ML-O-S	SC-2ML-Y-S	SC-2ML-G-S	SC-2ML-NA-S



Blue



Purple



White



Amber

Non-Sterile	SC-2ML-B	SC-2ML-PU	SC-2ML-W	SC-2ML-AM
Sterile	SC-2ML-B-S	SC-2ML-PU-S	SC-2ML-W-S	SC-2ML-AM-S

Screw Caps and Tubes

2ML

For the collection of clinical specimens containing virus, germs And for holding reagent

A variety colors of cap facilitate easy branding or organization.



Capture Screw Caps and Tubes

Monostick / Distick / Multistick

This product has a function that the cap part of the product can hold swab(s).

- Monostick Capture Caps can only catch one swab.
- Distick Capture Caps can catch two swabs.
- Multistick Capture Caps can catch up to four swabs.



Screw Caps and Tubes 5ML

For the collection of clinical specimens containing virus, germs And for holding reagent

5ML Screw Caps and Tubes, Non-Sterile

Cat. No.	Description	Cap Colors	Capture Type	Dimensions Ø x h	Packaging
SCT-5ML-RC	 Screw Caps and Tubes Non-Sterile	Red	Non	16.5mm x 82.5mm	2000 Tubes & Caps (500 Tubes x 4pk & 1000 Caps x 2pk)
SCTMO-5ML-RC		Red	Monostick		
SCTDI-5ML-RC		Red	Disitick		
SCTMU-5ML-RC		Red	Multistick		
SCT-5ML-WC		White	Non		
SCTMO-5ML-WC		White	Monostick		
SCTDI-5ML-WC		White	Disitick		
SCTMU-5ML-WC		White	Multistick		

5ML Screw Caps and Tubes, Sterile

Cat. No.	Description	Cap Colors	Capture Type	Dimensions Ø x h	Packaging
SCT-5ML-RC-S	 Screw Caps and Tubes Sterile by E-beam Irradiation	Red	Non	16.5mm x 82.5mm	2000 Tubes & Caps (500 Tubes x 4pk & 1000 Caps x 2pk)
SCTMO-5ML-RC-S		Red	Monostick		
SCTDI-5ML-RC-S		Red	Disitick		
SCTMU-5ML-RC-S		Red	Multistick		
SCT-5ML-WC-S		White	Non		
SCTMO-5ML-WC-S		White	Monostick		
SCTDI-5ML-WC-S		White	Disitick		
SCTMU-5ML-WC-S		White	Multistick		

Screw Caps and Tubes 12ML

For the collection of clinical specimens containing virus, germs And for holding reagent

12ML Screw Caps and Tubes, Sterile

Cat. No.	Description	Cap Colors	Capture Type	Dimensions Ø x h	Packaging
SCT-12ML-RC	Screw Caps and Tubes Non-Sterile	Red	Non	20mm x 103mm	1000 Caps and Tubes (1000 Caps/pk, 1pk/cs & 250 Tubes/pk, 4pk/cs)
SCT-12ML-GC		Green	Non		
SCT-12ML-WC		White	Non		



12ML Screw Caps and Tubes, Sterile

Cat. No.	Description	Cap Colors	Capture Type	Dimensions Ø x h	Packaging
SCT-12ML-RC-S	Screw Caps and Tubes Sterile by E-beam Irradiation	Red	Non	20mm x 103mm	1000 Caps and Tubes (1000 Caps/pk, 1pk/cs & 250 Tubes/pk, 4pk/cs)
SCT-12ML-GC-S		Green	Non		
SCT-12ML-WC-S		White	Non		



Screw Caps and Tubes 5ML

For the collection of clinical specimens containing virus, germs And for holding reagent

5ML Screw Caps, Non-Sterile

Cat. No.	Description	Cap Colors	Capture Type	Dimensions Ø x h	Packaging
SC-5ML-R	Screw Caps Non-Sterile 	Red	Non	16.5mm x 14.5mm	12000 Caps (2000 Caps/pk, 6 pks/cs)
SCMO-5ML-R		Red	Monostick		
SCDI-5ML-R		Red	Distick		
SCMU-5ML-R		Red	Multistick		
SC-5ML-W		White	Non		
SCMO-5ML-W		White	Monostick		
SCDI-5ML-W		White	Distick		
SCMU-5ML-W		White	Multistick		

5ML Screw Caps, Sterile

Cat. No.	Description	Cap Colors	Capture Type	Dimensions Ø x h	Packaging
SCT-5ML-RC-S	Screw Caps Sterile by E-beam Irradiation 	Red	Non	16.5mm x 14.5mm	12000 Caps (2000 Caps/pk, 6 pks/cs)
SCTMO-5ML-RC-S		Red	Monostick		
SCTDI-5ML-RC-S		Red	Distick		
SCTMU-5ML-RC-S		Red	Multistick		
SCT-5ML-WC-S		White	Non		
SCTMO-5ML-WC-S		White	Monostick		
SCTDI-5ML-WC-S		White	Distick		
SCTMU-5ML-WC-S		White	Multistick		

Screw Caps and Tubes

12ML

For the collection of clinical specimens containing virus, germs And for holding reagent

12ML Screw Caps, Non-Sterile

Cat. No.	Description	Cap Colors	Capture Type	Dimensions Ø x h	Packaging
SC-12ML-R	Screw Caps Non-Sterile	Red	Non	20mm x 14mm	6000 Caps (1000 Caps/pk, 6pk/cs)
SC-12ML-G		Green	Non		
SC-12ML-W		White	Non		

12ML Screw Caps, Sterile

Cat. No.	Description	Cap Colors	Capture Type	Dimensions Ø x h	Packaging
SC-12ML-R-S	Screw Caps Sterile by E-beam Irradiation	Red	Non	20mm x 14mm	6000 Caps (1000 Caps/pk, 6pk/cs)
SC-12ML-G-S		Green	Non		
SC-12ML-W-S		White	Non		



5ML, 12ML Screw Tubes

Cat. No.	Description	Colors	Dimensions Ø x h	Packaging
ST-5ML-C	Conical bottom, Self Standing, Non-Sterile	Clear	14mm x 81mm	500 tubes/pk, 4pk/cs
ST-5ML-C-S	Conical bottom, Self Standing, Sterile	Clear	14mm x 81mm	500 tubes/pk, 4pk/cs
ST-12ML-C	Conical bottom, Self Standing, Non-Sterile	Clear	17mm x 102mm	250 tubes/pk, 4pk/cs
ST-12ML-C-S	Conical bottom, Self Standing, Sterile	Clear	17mm x 102mm	250 tubes/pk, 4pk/cs





Bion Lifescience, Inc.

10F, 282, Sunhwagung-ro, Namyangju-si, Gyeonggi-do, Republic of Korea
bion@bionlifescience.com
Tel : +82-31-572-2272 Fax : +82-31-572-2273
www.bionlifescience.com